


*For Perforator Veins*

*Treatment of full spectrum venous insufficiency includes...*

# VNUS® ClosureRFS™ STYLET

*for Endovenous Radiofrequency Ablation of Incompetent Perforator Veins*



*Treating both superficial and perforator vein insufficiency results in*

- ***faster ulcer heal time***
- ***lower ulcer recurrence***

*than compression therapy alone.<sup>1,2</sup>*

## ***Product Features***

---

### ***> Efficacy***

- 84% to 93% vein occlusion<sup>3,4</sup>

### ***> Reliability***

- The only endovenous ablation device cleared by the FDA for the treatment of incompetent perforator veins
- Controlled, consistent thermal energy delivered to the vein wall

### ***> Quick Patient Recovery***

- Outpatient procedure
- Quick return to normal activities – usually within one day




**VNUS®**  
MEDICAL TECHNOLOGIES, INC.

# VNUS® ClosureRFS™

- > Unique design allows for various access methods under ultrasound guidance
- > Can be routinely used as a primary procedure, or adjunctive with the ClosureFAST™ catheter which treats superficial venous insufficiency



The ClosureRFS stylet is designed for use with the VNUS RFGPlus™ radiofrequency generators



VNUS ClosureRFS stylet treating perforator vein

#### REFERENCES:

- <sup>1</sup> Jamie R Barwell, Colin E Davies, Comparison of surgery and compression with compression alone in chronic venous ulceration (ESCHAR study): randomized controlled trial, THE LANCET, Vol 363, June 04
- <sup>2</sup> Nelzen O. Fransson I. True long-term healing and recurrence of venous leg ulcers following SEPS combined with superficial venous surgery: a prospective study. Eur J Vasc Endovasc Surg 34, 605-612 (2007)
- <sup>3</sup> Elias S, Peden E. Ultrasound-Guided Percutaneous Ablation for the Treatment of Perforating Vein Incompetence. Vascular 2007, Vol. 15, No. 5. pp.281-9.
- <sup>4</sup> RFS Registry: Data on file



ClosureRFS procedure pack also available

## VNUS ClosureRFS Ordering Information

Order#	Description	Device Size	Device Length
RFS2-6-12	RFS Stylet	6F	12cm
RFSP	Closure RFS Procedure Pack		

Compatible with .035" guidewire (with sharp removed)

> To order, call toll free: 1-888-797-VEIN (8346)



**VNUS**  
MEDICAL TECHNOLOGIES, INC.



**Corporate Office**  
VNUS Medical Technologies, Inc.  
5799 Fontanoso Way  
San Jose, CA 95138  
Tel: +1-408-360-7200  
Fax: +1-408-365-8480  
www.vnus.com



**Authorized Representative**  
VNUS Medical Technologies GmbH  
Wilhelm-Pfitzer-Str. 28  
70736 Fellbach  
Germany  
Tel: +49-(0)711-794809-0  
Fax: +49-(0)711-794809-3099

